

# Calibration of Black Carbon Instruments with the CPMA-Electrometer Reference Mass Standard (CERMS)

Joel C. Corbin, Gregory J. Smallwood, Prem Lobo  
NRC Metrology

Jason Olfert  
University of Alberta

24<sup>th</sup> ETH-Conference on Combustion Generated Nanoparticles  
22-24 June 2021 Virtual Meeting

# Acknowledgements

## National Research Council Canada

- Dan Clavel, Brett Smith, Simon-Alexandre Lussier



## Support from Transport Canada

- Wendy Bailey



Transport  
Canada

Transports  
Canada

## Rolls-Royce

- Mark Johnson and Matthew Parker



Rolls-Royce



# Outline

## Background

- techniques for measuring black carbon (BC)
- standards for aerosol instrument calibration
- the CPMA-electrometer reference mass standard (CERMS)

## Results

- calibration of a number of instruments
- calibration of an instrument with a number of sources

## Mass closure

- demonstration of mass closure when compared to elemental carbon (EC) from thermal optical analysis (TOA)

## Summary

# Techniques for measuring black carbon (BC)

## A non-exhaustive list of methods/instruments for measuring BC

Aerosol absorption methods

- CAPS PM<sub>ss</sub>a

Photoacoustic methods

- PAX, MSS, PASS

Filter-based methods

- MAAP, PSAP, Aethalometer

Evolved carbon methods

- TOA

Laser-induced incandescence

- LII 300, SP2

equivalent black carbon (eBC)


elemental carbon (EC)

refractory black carbon (rBC)

The reported values depend both on the characteristics of the source particles used in the calibration and the reference method applied, neither of which are unique or traceable

# Calibration of Condensation Particle Counters


$$N = I/Qe - N_{\text{multiply-charged}}$$



- charge
- classify
- compare to reference

# CPMA-electrometer reference mass standard (CERMS)

$$m_{total} = m_{+0} + M_1 I / Qe$$



- charge
- classify
- compare to reference

# Centrifugal Particle Mass Analyzer (CPMA)


## CPMA mass/charge setpoint depends on three traceable parameters:

- rotational speed (traceable)
- voltage (traceable)
- radial dimensions of the classifier
  - (optional traceable gauging certificate available from Cambustion)

Unlike a DMA, the setpoint does not depend upon the gas properties (such as viscosity and mean free path) or ambient conditions (temperature and pressure)

Classification not affected by particle morphology or composition

$$M_1 = \frac{M}{N_q} = \frac{eV}{\left(\frac{r_o + r_i}{2}\right)^2 \omega^2 \ln\left(\frac{r_o}{r_i}\right)}$$



# Benefits of CERMS

All aspects of the CERMS technique are traceable to the SI

Measurement of a well-defined quantity: total post-CPMA suspended PM mass

CERMS has low uncertainty<sup>1</sup>

- Repeatability was 1.1% ( $k = 1$ )
- Intermediate precision was 2.1% ( $k = 1$ )
- Resulting **uncertainty (coverage factor,  $k = 2$ )** was 3.0%


Rapid (<2 hours experiment time) with potential for same day determination of calibration factor

Additional benefits for determining aviation engine emissions BC mass

- Demonstrated to low concentrations ( $1 \mu\text{g}/\text{m}^3$ )
- on-site calibration
- gas turbine engine as the source of BC emissions



# CERMS apparatus for calibration with lab soot generator




# Calibration of LII 300 determined from CERMS with lab soot from MISG

## When referenced to denuded soot (nvPM):

$$\text{rBC} = (0.78 \pm 0.05) \cdot \text{nvPM}_{\text{CERMS}} + (0.2 \pm 0.1)$$

- rBC is refractory black carbon (LII 300)
- residuals are in the range of  $\pm 0.1 \mu\text{g}/\text{m}^3$
- fit intercept nearly zero (below the LOD of the instrument)
- note that for this source, the effect of the catalytic stripper (VPR) is negligible




# Calibration of PAX determined from CERMS with lab soot from MISG



## When referenced to nvPM,

$$\text{eBC} = (0.999 \pm 0.001) \cdot \text{nvPM}_{\text{CERMS}} - (1.05 \pm 0.02)$$


- eBC is equivalent black carbon (*i.e.*, calculated from absorption measurements with light-absorption-per-gram assumed as  $7.5 \text{ m}^2 \text{ g}^{-1}$  at  $550 \text{ nm}$ )
- $> 2 \mu\text{g eBC m}^{-3}$ , the calibration verifies the traditional optical PAX calibration
- $< 2 \mu\text{g eBC m}^{-3}$ , RH-related baseline problems become evident: CERMS provides the LOD of the instrument



# Calibration of LII 300 and MSS determined from CERMS with soot from large inverted burner




# CERMS Apparatus for calibration of a mass instrument on a gas turbine engine



# Calibration of MSS+ with CERMS using a gas turbine engine as the source – nvPM

## Calibration Factor


- Calibrated using the CERMS approach in under 2 hours, on-site at an OEM test cell
- For this MSS+, using the Gnome gas turbine engine as the source, the slope relating the LII response to nvPM is  $0.96 \pm 0.04$  with a small offset due to baseline effects



# Calibration of LII 300 with CERMS using a gas turbine engine as the source – nvPM

## Calibration Factor


- Calibrated using the CERMS approach in under 2 hours, on-site at an OEM test cell
- For this LII 300, using the Gnome gas turbine engine as the source, the slope relating the LII response to nvPM is  $0.63 \pm 0.02$



# Calibration of LII 300 with TOA using Gnome as the source – EC

## Calibration Factor

- Determined using the SAE ARP6320<sup>1</sup> method using EC collected on filters and thermal-optical analysis
- For the same LII 300, using the Gnome gas turbine engine as the source, the slope relating the LII response to EC is  $0.68 \pm 0.02$
- The difference in slopes reflects the difference in mass between nvPM and EC






# Demonstration of Closure Between CERMS and PM Mass from TOA

## Compare total PM mass from TOA to that from CERMS

- thermal-optical analysis (TOA) measures total carbon mass (TC), and operationally defines OC and EC fractions
- VOC contamination of the TOA filter was subtracted
- division of TC into OC and EC is often poorly defined for combustion samples
- as TC does not include all elements present in nvPM samples, it can be converted to total PM mass concentration
- %C is the mass fraction of carbon in the PM measured by elemental analysis (EA),  $f_C$



J. C. Corbin *et al.* (2020) *Closure between particulate mass measured ex situ by thermal-optical analysis and in situ by the CPMA-electrometer reference mass system*, *Aerosol Science and Technology*, 54, 1293-1309.

doi:[10.1080/02786826.2020.1788710](https://doi.org/10.1080/02786826.2020.1788710)

# Summary

## **CERMS (CPMA-Electrometer Reference Mass Standard) has been demonstrated for rapid, low uncertainty calibration of a range of BC instruments using a number of BC sources**

This study has successfully demonstrated operation of a CERMS calibration of mass instruments using a gas turbine engine as the source

- Rapid (<2 hours experiment time) with same-day determination of calibration factor
- May be performed on-site
- Sensitive to very low concentrations, relevant to nvPM mass concentrations from modern gas turbine engines
- Instrument responses are linear with respect to CERMS regardless of instrument detection methodology or source

# THANK YOU

**Dr. Gregory J. Smallwood**

**Principal Research Officer**

**NRC Metrology**

phone: 613-993-1391

e-mail: [greg.smallwood@nrc-cnrc.gc.ca](mailto:greg.smallwood@nrc-cnrc.gc.ca)


# The CPMA



The **CPMA** classifies particles by mass/charge ratio

An aerosol **electrometer** measures total particulate charge,

so a **CPMA**-electrometer measures total particulate mass



# What about uncharged particles?

They are generally close to zero, even so, they can be directly measured using an electrostatic precipitator to remove charged particles

Particle source



Unipolar charger



Centrifugal Particle Mass Analyzer



Aerosol electrometer



Measures the current,  $I$ , which is proportional to number of charges detected

pump

Electrostatic Precipitator



Mass instrument to calibrate

